

21<sup>ST</sup>



# PGDIS CONFERENCE



6-8 May 2024  
Kuala Lumpur  
Malaysia

**PGT and  
BEYOND...**




# Sample swap and contamination detection in PGT-A

Pedro Echave, PhD





# The problem of contamination in PGT-A



PGT-A is a widely used technique during in vitro fertilization procedures, serving as a decision-making tool. Detecting contamination during PGT-A is crucial to prevent misclassification of the embryos.

## Implications\*

- Discard potentially viable embryos
- Transfer of aneuploid embryos wrongly classified as mosaic, potentially leading to increased rates of implantation failure and miscarriage.

\*Clark, G., et al. (2023). *Human Reproduction*, 38(Supplement\_1), dead093.337

# Lab practices to mitigate contamination risk

- 1. Negative Controls:** Implement always negative controls during sample processing. Limited by the fact that cannot evaluate contamination in the actual tube containing the biopsy.
- 2. Laboratory Practices:** Adhere to strict laboratory practices to minimize the risk of contamination. This includes maintaining a clean environment, proper gowning procedures, and sterilization of equipment.
- 3. Barrier Technology:** Utilize barrier technology to prevent cross-contamination. This can involve physical barriers such as laminar flow hoods or biological safety cabinets.
- 4. Personnel Training:** Train laboratory personnel thoroughly on contamination control measures. Regular training and awareness programs can help maintain a high level of vigilance.

## Extent of the problem\*


On average 0.4% variable from clinic to clinic

\*Clark, G., et al. (2023). *Human Reproduction*, 38(Supplement\_1), dead093.337



# Contamination detection with WGA-based method

## PG-Seq™ Rapid v2




- Ready-mixed reagents
- Minimal tube opening and pipetting steps
- No normalization required

### Degenerate Oligo PCR-based method

- Amplifies nuclear gDNA
- **Amplifies mtDNA**


# mtDNA amplification



Real D5 embryo sample processed with PG-Seq<sup>TM</sup> Rapid v2

Alignment visualized with IGV to illustrate coverage

- FastQ generated in the standard PG-Seq Rapid<sup>TM</sup> v2 workflow can be used to obtain information about mtDNA sequence of the sample
- Percentage of reads aligning to mitochondrial genome range from 0.4% to 20.2%, with average of 4.8% (n = 1140)





## Pros

### 1. Unique Inheritance Pattern:

- Unlike nuclear DNA, mtDNA is **maternally inherited**.
- Each embryo inherits its mtDNA exclusively from the mother.
- By analyzing mtDNA variants, we can trace the maternal lineage and uniquely identify embryos.

### 2. Contamination Detection:

- External contamination (e.g., from lab personnel or sample handling) can introduce foreign mtDNA.
- By comparing the embryo's mtDNA profile with expected maternal mtDNA, we can identify any discrepancies.
- Contaminated samples can be flagged, ensuring accurate results.


### 3. Confirmation of Sibling Relationships:

- When multiple embryos are tested, mtDNA confirms whether they share the same maternal lineage.
- It ensures that tested embryos are indeed full genetic siblings.

## Cons

### 1. Maternal DNA contamination

- Maternal and embryonic mtDNA are identical,
- If contamination from cumulus cells is present in the biopsy, it won't be detected by this method




[http://www.ijmr.org.in/viewimage.asp?img=IndianJMedRes\\_2017\\_146\\_3\\_341\\_223630\\_f2.jpg](http://www.ijmr.org.in/viewimage.asp?img=IndianJMedRes_2017_146_3_341_223630_f2.jpg)

# Proof of principle

- Custom bioinformatic analysis was developed that extracts mtDNA SNV then compares similarity between samples.
- Known pathogenic variants and low depth SNV (depth <2) excluded.
- 1140 samples were processed and sequenced, with 5.4M total read on average.
- After downsampling to 500,000 reads, fastq were re-run through the pipeline

	5.4M avg read	0.5M read
Correctly grouped	1137 (99.73%)	1128 (98.94%)
Incorrectly grouped together	1	5
Incorrectly grouped separately	0	2
Not enough mtDNA reads for analysis	2	5








- Cloud-based software is free and to be used as an additional tool along with the PG-find software\*.
- Output is an excel containing sample groupings
- Can be used to detect exogenous contamination and sample swaps\*\*.

Group	Samples															
G1	01															
G2	02															
G3	03	04	05	06	07	08	09	10								
G4	11	12	13													
G5	14															
G6	15	16	17													
G7	18	19	20													
G8	21															
G9	22	23	24													
G10	25															
G11	26	27	28													
G12	29	30	31	32	33	34	35	36	37	38	39					
G13	40															
G14	41															
G15	42	43														
G16	44	45														
G17	46	47	48	49	50											

\*Server currently located in USA

\*\*Percentage of sample swapping taking place in current PGT-A labs is unknown to the best of our knowledge, but consequence is same as contamination: misclassification

# mtDetect™ Web App





# THANK YOU

**PGT and  
BEYOND...**